

Product Datasheet

Anti-Shiga toxin (Stx-2b) V7587

Product Name	Anti-Shiga toxin (Stx-2b)
Catalog Number	V7587
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Shiga toxin type II - B subunit (Stx- 2b)
Shelf Life	Four years from date of manufacture