

Product Datasheet

Anti-Cryptosporidium parvum V7631

Product Name	Anti-Cryptosporidium parvum
Catalog Number	V7631
Tested Applications	ELISA, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with intact oocysts of <i>C. parvum</i>
Specificity	Intact oocysts of <i>C. parvum</i>
Shelf Life	Four years from date of manufacture