

Product Datasheet

Anti-Cholera toxin (Beta subunit) V7954

Product Name	Anti-Cholera toxin (Beta subunit)
Catalog Number	V7954
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with intact toxin. Non-reactive with E. Coli heat labile toxin
Specificity	B subunit of cholera toxin
Matched Pairs	ELISA pairs (capture-conjugate) : V7954-V7953, V7954-V7955
Shelf Life	Four years from date of manufacture