

Product Datasheet

Anti-Vaccinia virus V8112

Product Name	Anti-Vaccinia virus
Catalog Number	V8112
Tested Applications	LF, ELISA, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with free virus, uninfected cells, Lister and NYCBH strains of vaccinia
Specificity	Free virus + infected cells
Matched Pairs	LF pairs (capture-conjugate) : V8112-V8114
Shelf Life	Four years from date of manufacture