

Product Datasheet

Anti-h ApoA1 2002 SPTN-5 100264

Product Name	Anti-h ApoA1 2002 SPTN-5
Catalog Number	100264
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	KA= 1 x 10 ⁹ l/M; KD= 1 x 10 ⁻⁹ M
Associated Products	Recombinant ApoA1, Medix Biochemica, 710012, 710033, 710032
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.7-7.7
Cross Reactivity	N/D
Specificity	Antibody recognizes human apolipoprotein A1 (ApoA1)
Shelf Life	36 months