

Product Datasheet

Anti-Human CD20 Monoclonal Antibody, Azide Free Clone B-H20 854.150

Product Name	Anti-Human CD20 Monoclonal Antibody, Azide Free Clone B-H20
Catalog Number	854.150
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Functional assay, FC
Immunogen	Chronic Lymphoid Leukaemia cells (CLL)
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the B1, Bp35 antigen, a 33-37 kDa protein
Biological Activity	Cytotoxic in association with complement on CD20 cells
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells