

Product Datasheet

Anti-Feline Leukemia virus (p27) V7228

Product Name	Anti-Feline Leukemia virus (p27)
Catalog Number	V7228
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	FeLv p27
Matched Pairs	Pairs for p27 detection (capture-conjugate) : V7228-V7218, V7228-V7228
Shelf Life	Four years from date of manufacture