

Product Datasheet

Anti-Canine Parvovirus V7351

Product Name	Anti-Canine Parvovirus
Catalog Number	V7351
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with CPV and Feline Panleukopenia virus
Specificity	CPV and Feline Panleukopenia virus
Matched Pairs	Pairs for antigen detection (capture-conjugate) : V7352-V7351, V7341-V7351, V7353-V7351
Shelf Life	Four years from date of manufacture