

## Product Datasheet

### Anti-Canine Parvovirus V7351

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Anti-Canine Parvovirus  |
| <b>Catalog Number</b>      | V7351   |
| <b>Tested Applications</b> | ELISA, IFA  |
| <b>Brand</b>               | ViroStat  |
| <b>Form/Appearance</b>     | Purified from mouse ascites by protein A chromatography                                 |
| <b>Concentration</b>       | ≥ 1 mg/ml   |
| <b>Storage</b>             | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.             |
| <b>Isotype</b>             | IgG2a   |
| <b>Clonality</b>           | Monoclonal  |
| <b>Purity</b>              | ≥ 90%   |
| <b>Buffer</b>              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide                   |
| <b>Cross Reactivity</b>    | Reactive with CPV and Feline Panleukopenia virus  |
| <b>Specificity</b>         | CPV and Feline Panleukopenia virus  |
| <b>Matched Pairs</b>       | Pairs for antigen detection (capture-conjugate) : V7352-V7351, V7341-V7351, V7353-V7351 |
| <b>Shelf Life</b>          | Four years from date of manufacture   |