

## Product Datasheet

### Anti-Canine Parvovirus V7354

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Anti-Canine Parvovirus  |
| <b>Catalog Number</b>      | V7354   |
| <b>Tested Applications</b> | ELISA, IFA  |
| <b>Brand</b>               | ViroStat  |
| <b>Form/Appearance</b>     | Purified from mouse ascites by protein A chromatography                     |
| <b>Concentration</b>       | ≥ 1 mg/ml   |
| <b>Storage</b>             | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Purity</b>              | ≥ 90%   |
| <b>Buffer</b>              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide       |
| <b>Cross Reactivity</b>    | Reactive with CPV. Non-reactive with Feline Panleukopenia virus             |
| <b>Specificity</b>         | CPV   |
| <b>Shelf Life</b>          | Four years from date of manufacture   |