

Product Datasheet

Hepatitis C virus (NS3) Antigen, Recombinant V8801

Product Name Hepatitis C virus (NS3) Antigen, Recombinant

Catalog Number V8801

Species/Host E. coli

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C. Avoid freeze-thaws.

Note No fusion partner present

Purity ~90% by SDS-PAGE

Buffer 25 mM Tris, pH 8.3 + 0.1% SDS

Specificity NS3 of HCV, 1230aa-1501aa

Molecular Weight 13 kDa