

## Product Datasheet

### Hepatitis C virus (core) Antigen, Recombinant V8802

**Product Name** Hepatitis C virus (core) Antigen, Recombinant

**Catalog Number** V8802

**Species/Host** E. coli

**Tested Applications** ELISA

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** ≤ -20°C

**Note** No fusion partner present

**Purity** ~95% by SDS-PAGE

**Buffer** 10 mM PBS, pH 7.4

**Specificity** HCV core

**Molecular Weight** 66 kDa