

Product Datasheet

Rubella virus (E1) Antigen, Recombinant V8806

Product Name Rubella virus (E1) Antigen, Recombinant

Catalog Number V8806

Species/Host E. coli

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note No fusion partner present

Purity ~95% by SDS-PAGE

Buffer 25 mM Tris-HCl, 0.1% SDS, pH 8.3

Specificity E1

Molecular Weight 53 kDa