

## Product Datasheet

### Rubella virus (E2) Antigen, Recombinant V8807

**Product Name** Rubella virus (E2) Antigen, Recombinant

**Catalog Number** V8807

**Species/Host** E. coli

**Tested Applications** ELISA

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** ≤ -20°C

**Note** No fusion partner present

**Purity** ~95% by SDS-PAGE

**Buffer** 25 mM Tris-HCl, 0.1% SDS, pH 8.3

**Specificity** E2

**Molecular Weight** 38 kDa