

Product Datasheet

Hepatitis C virus (NS3) Antigen, Recombinant V8819

Product Name	Hepatitis C virus (NS3) Antigen, Recombinant
Catalog Number	V8819
Species/Host	E. coli
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	GST
Purity	~90% by SDS-PAGE
Buffer	10 mM Tris, pH 8.5 + 10 mM NaCl + 0.5M EDTA + inhibitors
Specificity	NS3 of HCV, 272 aa