

Product Datasheet

Toxoplasma gondii (SAG1) Antigen, Recombinant V8821

Product Name	Toxoplasma gondii (SAG1) Antigen, Recombinant
Catalog Number	V8821
Species/Host	E. coli
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	6X-His
Purity	~95% by SDS-PAGE
Buffer	25 mM Tris-HCl, pH 8.3 + 1.5M Urea + 50% glycerol
Specificity	SAG1 (p30), aa45-aa198