

Product Datasheet

Anti-Human CD20 Monoclonal Antibody, Unconjugated Clone B-H20 854.153.020

Product Name	Anti-Human CD20 Monoclonal Antibody, Unconjugated Clone B-H20
Catalog Number	854.153.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Chronic Lymphoid Leukaemia cells (CLL)
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the B1, Bp35 antigen, a 33-37 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells