

Product Datasheet

Hepatitis C virus (core+NS3+NS4+NS5) Antigen, Recombinant V8823

Product Name Hepatitis C virus (core+NS3+NS4+NS5) Antigen, Recombinant

Catalog Number V8823

Species/Host E. coli

Tested Applications LF

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note No fusion partner present

Purity ~95% by SDS-PAGE

Buffer PBS, 7.2 + 125mM arginine

Specificity HCV 4th generation, core + NS3 + NS4 + NS5

Molecular Weight 20 kDa