

Product Datasheet

Dengue virus (Env, Serotypes 1-4) Antigen, Recombinant V8824

Product Name Dengue virus (Env, Serotypes 1-4) Antigen, Recombinant

Catalog Number V8824

Species/Host E. coli

Tested Applications LF, ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note 6X-His

Purity ~95% by SDS-PAGE

Buffer PBS, pH 7.4

Specificity Envelope antigen, contains serotype common regions

Molecular Weight 19 kDa