

## Product Datasheet

### Ebola virus (NP) Antigen, Recombinant V8828

<b>Product Name</b>	Ebola virus (NP) Antigen, Recombinant
<b>Catalog Number</b>	V8828
<b>Species/Host</b>	Pichia pastoris
<b>Tested Applications</b>	LF, ELISA, WB
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	≤ -20°C
<b>Note</b>	No fusion partner present
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	PBS, 7.2 + 0.02% sodium azide
<b>Specificity</b>	Nucleoprotein (NP), C terminus, Sudan strain
<b>Molecular Weight</b>	16 kDa