

## Product Datasheet

### Hepatitis C virus (core) Antigen, Recombinant V8903

<b>Product Name</b>	Hepatitis C virus (core) Antigen, Recombinant
<b>Catalog Number</b>	V8903
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	LF, ELISA, WB
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	≤ -20°C
<b>Note</b>	β-galactosidase (114 kDa)
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	0.05M Tris-saline, pH 8.0 + 6M urea
<b>Specificity</b>	Core antigen of HCV, truncated from C-terminus
<b>Molecular Weight</b>	133 kDa