

## Product Datasheet

### Hepatitis C virus (NS4) Antigen, Recombinant V8904

<b>Product Name</b>	Hepatitis C virus (NS4) Antigen, Recombinant
<b>Catalog Number</b>	V8904
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	ELISA, WB
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	≤ -20°C
<b>Note</b>	β-galactosidase (114 kDa)
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	8 M Urea, 20 mM Tris-HCl pH 8.0
<b>Specificity</b>	NS4 of HCV, NS4a+b; a.a. 1658 thru a.a. 1863 (19 kDa)
<b>Molecular Weight</b>	146 kDa