

## Product Datasheet

### Human Immunodeficiency Virus 1 (gp41) Antigen, Recombinant V8908

**Product Name** Human Immunodeficiency Virus 1 (gp41) Antigen, Recombinant

**Catalog Number** V8908

**Species/Host** Pichia pastoris

**Tested Applications** ELISA

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** 2-8°C

**Note** No fusion partner present

**Purity** ~95%

**Buffer** 10 mM Na<sub>2</sub>CO<sub>3</sub>/10 mM EDTA/14 mM β-ME/0.05% Tween 20

**Specificity** Envelope, gp41

**Molecular Weight** 26 kDa