

Product Datasheet

Human Immunodeficiency Virus 1 (p24) Antigen, Recombinant V8909

Product Name Human Immunodeficiency Virus 1 (p24) Antigen, Recombinant

Catalog Number V8909

Species/Host E. coli

Tested Applications ELISA, WB

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note β-galactosidase (114 kDa)

Purity ~95%

Buffer 0.02 M Sodium phosphate, pH 7.2

Specificity Gag, p24

Molecular Weight 148 kDa