

Product Datasheet

Human Immunodeficiency Virus 2 (gp36) Antigen, Recombinant V8911

Product Name Human Immunodeficiency Virus 2 (gp36) Antigen, Recombinant

Catalog Number V8911

Species/Host E. coli

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note No fusion partner present

Purity ~95%

Buffer 10 mM Na₂CO₃/10 mM EDTA/14 mM β-ME/0.05% Tween 20

Specificity Envelope, gp36

Molecular Weight 13.5 kDa