

Product Datasheet

Hepatitis B virus (pre-S1) Antigen, Recombinant V8913

Product Name	Hepatitis B virus (pre-S1) Antigen, Recombinant
Catalog Number	V8913
Species/Host	E. coli
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	No fusion partner present
Purity	~95% by SDS-PAGE
Buffer	PBS, pH 7.2
Specificity	Pre-S1 antigen, a.a. #1-a.a. #108 of adw pre-S region