

## Product Datasheet

### Hepatitis B virus Core (HBcAg) Antigen, Recombinant V8914

**Product Name** Hepatitis B virus Core (HBcAg) Antigen, Recombinant

**Catalog Number** V8914

**Species/Host** E. coli

**Tested Applications** ELISA, WB

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** ≤ -20°C

**Note** No fusion partner present

**Purity** ~95% by SDS-PAGE

**Buffer** 50 mM Tris, 1mM EDTA, pH 7.5 + 150 mM NaCl

**Specificity** Core antigen, 180 a.a.

**Molecular Weight** 17 kDa