

Product Datasheet

Hepatitis B virus (HBeAg) Antigen, Recombinant V8915

Product Name	Hepatitis B virus (HBeAg) Antigen, Recombinant
Catalog Number	V8915
Species/Host	Pichia pastoris
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	No fusion partner present; excellent correlation with a commercial EBNA ELISA when used as a coating antigen @ 10 ug/ml
Purity	~95% by SDS-PAGE
Buffer	50 mM Sodium phosphate, pH 6.8 + 5% glycerol
Specificity	"e" antigen, 160 a.a.
Molecular Weight	35 kDa