

## Product Datasheet

### Epstein-Barr virus (EBNA-1) Antigen, Recombinant V8916

<b>Product Name</b>	Epstein-Barr virus (EBNA-1) Antigen, Recombinant
<b>Catalog Number</b>	V8916
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	2-4°C
<b>Note</b>	β-galactosidase (114 kDa)
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	Phosphate buffered saline, pH 7.2 + 0.1% sodium azide
<b>Specificity</b>	EBNA-1, truncated to 340 aa, glycine-alanine repeats
<b>Molecular Weight</b>	135 kDa