

## Product Datasheet

### Hepatitis C virus (core) Antigen, Recombinant V8924

<b>Product Name</b>	Hepatitis C virus (core) Antigen, Recombinant
<b>Catalog Number</b>	V8924
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	2-8°C
<b>Note</b>	β-galactosidase (114 kDa)
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	20 mM Tris-HCl, pH 8 + 10 mM β-ME + 8M urea
<b>Specificity</b>	HCV core, a.a. 2 thru a.a. 192 (22 kDa)
<b>Molecular Weight</b>	18 kDa