

## Product Datasheet

### Hepatitis B virus Core (HBcAg) Antigen, Recombinant V8927

**Product Name** Hepatitis B virus Core (HBcAg) Antigen, Recombinant

**Catalog Number** V8927

**Species/Host** E. coli

**Tested Applications** ELISA

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** ≤ -20°C

**Note** No fusion partner present

**Purity** ~95% by SDS-PAGE

**Buffer** PBS/glycerol (50%)

**Specificity** Core antigen of Hepatitis B virus, a.a. 1-183

**Molecular Weight** 18 kDa