

Product Datasheet

Hepatitis B virus Core (HBcAg) Antigen, Recombinant V8929

Product Name Hepatitis B virus Core (HBcAg) Antigen, Recombinant

Catalog Number V8929

Species/Host *S. cerevisiae*

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note No fusion partner present

Purity ~95% by SDS-PAGE

Buffer PBS/glycerol (50%)

Specificity Core antigen of Hepatitis B virus, a.a. 1-183

Molecular Weight 24 kDa