

Product Datasheet

Hepatitis B virus (HBsAg) Antigen, Recombinant V8930

Product Name Hepatitis B virus (HBsAg) Antigen, Recombinant

Catalog Number V8930

Species/Host S. cerevisiae

Tested Applications ELISA, WB

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage 2-8°C

Note No fusion partner present

Purity ~95% by SDS-PAGE

Buffer PBS, pH 7.2

Specificity Surface antigen, ay subtype

Molecular Weight 24 kDa