

## Product Datasheet

### Hepatitis B virus (HBsAg) Antigen, Recombinant V8931

**Product Name** Hepatitis B virus (HBsAg) Antigen, Recombinant

**Catalog Number** V8931

**Species/Host** Pichia pastoris

**Tested Applications** ELISA, WB

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** 2-8°C

**Note** No fusion partner present

**Purity** ~95% by SDS-PAGE

**Buffer** PBS, pH 7.2

**Specificity** Surface antigen, ay subtype

**Molecular Weight** 22-28 kDa; multiple glycosylated bands