

Product Datasheet

Hepatitis E virus (ORF 2) Antigen, Recombinant V8942

Product Name Hepatitis E virus (ORF 2) Antigen, Recombinant

Catalog Number V8942

Species/Host E. coli

Tested Applications ELISA, WB

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note β-galactosidase

Purity ~95% by SDS-PAGE

Buffer 20 mM Tris, pH 8.0 + 8M urea + 10 mM BME

Specificity ORF2 antigen of HEV, a.a.633-a.a.659