

## Product Datasheet

### Hepatitis E virus (ORF 2) Antigen, Recombinant V8943

**Product Name** Hepatitis E virus (ORF 2) Antigen, Recombinant

**Catalog Number** V8943

**Species/Host** E. coli

**Tested Applications** ELISA

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** ≤ -20°C

**Note** β-galactosidase

**Purity** ~95% by SDS-PAGE

**Buffer** 20 mM Tris, pH 8.0 + 8M urea + 10 mM BME

**Specificity** ORF2 antigen of HEV, a.a.403- a.a.461