

Product Datasheet

Treponema pallidum (p17) Antigen, Recombinant V8946

Product Name Treponema pallidum (p17) Antigen, Recombinant

Catalog Number V8946

Species/Host E. coli

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note β-galactosidase

Purity ~95% by SDS-PAGE

Buffer 20 mM Tris, pH 8.0 + 8M urea + 10 mM BME

Specificity p17, full length

Molecular Weight ~ 160 kDa