

Product Datasheet

Treponema pallidum (TnpA) Antigen, Recombinant V8948

Product Name Treponema pallidum (TnpA) Antigen, Recombinant

Catalog Number V8948

Species/Host Pichia pastoris

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note No fusion partner present

Purity ~95% by SDS-PAGE

Buffer 20 mM Tris, pH 8.0 + 8M urea + 10 mM BME

Specificity TnpA, full length

Molecular Weight 14 kDa