

Product Datasheet

Human Immunodeficiency Virus 2 (gp36) Antigen, Recombinant V8949

Product Name Human Immunodeficiency Virus 2 (gp36) Antigen, Recombinant

Catalog Number V8949

Species/Host *S. cerevisiae*

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note Human superoxide dismutase

Purity ~95% by SDS-PAGE

Buffer 50 mM Tris, pH 8.0 + 0.5M NaCl + 25 mM imidazole + 6M urea

Specificity Envelope, portion of gp36

Molecular Weight 35 kDa