

Product Datasheet

Herpes Simplex Virus 1 (gG) Antigen, Recombinant V8950

Product Name Herpes Simplex Virus 1 (gG) Antigen, Recombinant

Catalog Number V8950

Species/Host S. cerevisiae

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note Human superoxide dismutase

Purity ~95% by SDS-PAGE

Buffer 20 mM Phosphate, pH 7.4 + 6M urea

Specificity Glycoprotein gG, a.a. 1 thru a.a. 90