

Product Datasheet

Herpes Simplex Virus 2 (gG) Antigen, Recombinant V8951

Product Name Herpes Simplex Virus 2 (gG) Antigen, Recombinant

Catalog Number V8951

Species/Host E. coli

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note Terminal 6-his

Purity ~95% by SDS-PAGE

Buffer 50 mM Phosphate, pH 7.0 + 50mM DTT

Specificity Glycoprotein gG, a.a. 343 thru a.a. 649

Molecular Weight 21 kDa