

Product Datasheet

Hepatitis B virus (Pre-S1 + S2) Antigen, Recombinant V8952

Product Name	Hepatitis B virus (Pre-S1 + S2) Antigen, Recombinant
Catalog Number	V8952
Species/Host	E. coli
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	Terminal 6-his
Purity	~95% by SDS-PAGE
Buffer	PBS, pH 7.2
Specificity	Pre-S1+S2 antigen, a.a. #1-a.a. #174 of adw pre-S1 region
Molecular Weight	7.1 kDa