

Product Datasheet

Feline Immunodeficiency Virus (core) Antigen, Recombinant V8963

Product Name Feline Immunodeficiency Virus (core) Antigen, Recombinant

Catalog Number V8963

Species/Host E. coli

Tested Applications LF, ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage 2-8°C

Note Maltose binding protein (truncated)

Purity ~90% by SDS-PAGE

Buffer Bicarbonate buffer, pH 8.5 + 0.02% azide

Specificity Core

Molecular Weight 37 kDa