

## Product Datasheet

### Plasmodium falciparum (protein II) Antigen, Recombinant V8964

**Product Name** Plasmodium falciparum (protein II) Antigen, Recombinant

**Catalog Number** V8964

**Species/Host** E. coli

**Tested Applications** ELISA

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** ≤ -20°C

**Note** GST

**Purity** ~90% by SDS-PAGE

**Buffer** 20 mM Carbonate, pH9.6 + 0.1% sodium azide

**Specificity** Histidine rich protein-2, full length

**Molecular Weight** 32 kDa