

Product Datasheet

Cytomegalovirus (gB) Antigen, Recombinant V8965

Product Name	Cytomegalovirus (gB) Antigen, Recombinant
Catalog Number	V8965
Species/Host	E. coli
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	No fusion partner present
Purity	~95% by SDS-PAGE
Buffer	25 mM Tris, pH 8.0 + 1 mM EDTA + 50% glycerol
Specificity	gB, aa.11 to aa.67
Molecular Weight	30 kDa