

Product Datasheet

Cytomegalovirus (p28) Antigen, Recombinant V8966

Product Name	Cytomegalovirus (p28) Antigen, Recombinant
Catalog Number	V8966
Species/Host	E. coli
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	Maltose binding protein (truncated)
Purity	~95% by SDS-PAGE
Buffer	25 mM Tris, pH 8.0 + 1 mM EDTA + 50% glycerol
Specificity	p28 (UL99), a.a. 130 to aa. 160
Molecular Weight	~25 kDa