

Product Datasheet

Plasmodium (aldolase) Antigen, Recombinant V8970

Product Name	Plasmodium (aldolase) Antigen, Recombinant
Catalog Number	V8970
Species/Host	E. coli
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	GST
Purity	~95% by SDS-PAGE
Buffer	10 mM Tris, pH 8.0 + 10 mM NaCl + 0.02% ProClin 300
Specificity	Aldolase, 368 a.a.
Molecular Weight	38 kDa