

Product Datasheet

Hepatitis C virus (NS5a) Antigen, Recombinant V8994

Product Name Hepatitis C virus (NS5a) Antigen, Recombinant

Catalog Number V8994

Species/Host E. coli

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage 2-8°C

Note GST

Purity ~95% by SDS-PAGE

Buffer 50 mM Tris, pH 8.0 + EDTA

Specificity NS5 antigen of HCV, NS5a region

Molecular Weight 68 kDa