

Product Datasheet

Human Immunodeficiency Virus 1 (Subtype B gp 41) Antigen, Recombinant V9000

Product Name	Human Immunodeficiency Virus 1 (Subtype B gp 41) Antigen, Recombinant
Catalog Number	V9000
Species/Host	E. coli
Tested Applications	LF, ELISA, WB
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	His tag
Purity	~95% by SDS-PAGE
Buffer	20 mM Na-carbonate, pH 9.6
Specificity	Ectodomain Subtype B, gp41
Molecular Weight	13 kDa