

Product Datasheet

Human Immunodeficiency Virus 1 (subtype C p24) Antigen, Recombinant V9002

Product Name	Human Immunodeficiency Virus 1 (subtype C p24) Antigen, Recombinant
Catalog Number	V9002
Species/Host	E. coli
Tested Applications	LF, ELISA, WB
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	His tag
Purity	~95% by SDS-PAGE
Buffer	PBS, pH 7.4
Specificity	Subtype C clone, p24
Molecular Weight	17.0 kDa